


**New**

# ALASTA<sup>®</sup> PRO

Thin Nitrile Exam Gloves



**Putting on gloves has never been easier!**



- *Gloves glide on damp hands easily*
- *Soft, form fitting blue nitrile*
- *Tested against fentanyl*
- *Value Packed 200 gloves/box*

**3 mil**



**SAY GOODBYE TO THE FRUSTRATION OF PUTTING ON GLOVES!**

Introducing...Glide Assist™ Technology. Thanks to our proprietary inner coating, your hands glide right in, no matter how damp they are. Get your gloves on right the first time - no wasted time, no wasted money.\*

**7/10**

*Respondents perceived that the glove was easier to put on and had less resistance.\*\**

\*\*Per independent testing

**Your Trusted Choice in Gloves - Protecting Hands and People, Made Easy**



800-523-2055 | DASHMEDICAL.COM

A new addition to our Alasta<sup>®</sup> family. Introducing Alasta<sup>®</sup>PRO with Glide Assist<sup>™</sup> technology! Putting on gloves has never been easier. Glide Assist<sup>™</sup> eliminates the hassle of putting on gloves, so you can get to work faster. Made with Alasta Soft Fit nitrile for all day comfort and strength for peace of mind.



#### PRODUCT DETAILS

Material	Nitrile
Interior	Powder Free
Grade	Medical Examination
Watertight Test	1.0 AQL
Color	Blue
Texture	Finger

#### PHYSICAL PROPERTIES (TYPICAL VALUES\*)

Palm Thickness	0.07 mm   2.8 mil
Finger Thickness	0.10 mm   3.9 mil
Length	240 mm   9.5 in
Tensile Strength (MPA)	17
Elongation %	500

#### ORDER INFORMATION

Product Number	ALGA200
Available Sizes	XS, S, M, L, XL, XXL
Glove Count	200/box 10 boxes/case

#### CERTIFICATIONS

Chemotherapy Drug Tested <sup>1</sup>	Yes
Tested Against Fentanyl	Yes
PFAS Free	Yes
BPA Free	Yes
FDA Food Contact	Yes

#### PRODUCT CONFORMANCE

ASTM D6319, ASTM D6978 (Chemotherapy & Fentanyl), ASTM D5151, ASTM D6124, FDA 21 CFR 177

#### RECOMMENDED INDUSTRIES

Corrections, Dental, EMS, Healthcare, Veterinarian, Beauty and Industrial Application

PERFORMANCE LEVEL	LIGHT DUTY	HIGH TACTILITY	MID DUTY	HEAVY DUTY	HIGH RISK
CUSTOMER NEED/ APPLICATION	Basic Barrier Coverage	Best Touch Sensitivity	Broad Range of Uses	Increased Durability	Longer Cuff Length
	Lowest cost option	Delivers a second skin feel	Perfect for a variety of applications	Great for more rugged applications	Used in higher risk settings

\*Glide Assist<sup>™</sup> is a safe water-based inner coating that reduces friction allowing gloves to easily slide over the user's hand.

See DASH<sup>®</sup> in Dental Product Shopper & Forbes



**DASH**  
MEDICAL GLOVES

800-523-2055 | DASHMEDICAL.COM

DASH\_SS\_ALGA200\_07.2025

Contact DASH Medical for a list of tested Chemotherapy drugs  
Contact DASH Medical for a copy of the independent testing